

SPECTROPHOTOMETER CM-3600d

Built for Precision, Priced for Economy



New Enhanced Performance Technology teams patented performance features with simplified design and advanced cost-reducing manufacturing techniques.

The result:

- A highly accurate, reliable, rugged spectrophotometer.
- Versatile enough for all colorimetric applications.
- Simplified operation.

And

■ LOW PRICE

Technology

KONICA MINOLTA Innovative Optical System "World first" technology:

- Numerical Gloss Control measurement within a few seconds.
- Numerical UV Control for calibration and measurements.
- Soft-flash mode to avoid triplet adsorption in fluorescent materials.

The Innovative CM-3600d Spectrophotometer Brings You Highest Quality Levels in Color Measurement along with Unsurpassed Versatility, User Comfort and Low Cost!

Enhanced Performance Technology Konica Minolta Innovative Optical System


The CM-3600d is equipped with Minolta Innovative Optical System technology. For users, this means high accuracy and repeatability, improved performance, additional features, simple operation, and - low, affordable pricing.


This technology also provides:

1.Numerical Gloss Control

Instead of using a mechanically driven gloss trap, the CM-3600d is equipped with a patented numerical Gloss (SCI/SCE) control system. By sequentially firing two flashes, within a few seconds the system provides both SCI and SCE values for each sample.


SCI ······ specular component included SCE···· specular component excluded






2.Numerical UV Control


To replace the time-consuming traditional UV measurement system that utilizes moving filters, the CM-3600d introduces yet another patented numerical calculation system, never seen before: The measured values of two sequentially fired flashes, one with full UV energy and one with UV cut-off filter at either 400 or 420 nm are combined to obtain the spectral characteristics, and the respective whiteness and tint value of any UV-activated fluorescence sample. To avoid triplet effect on FWA treated samples, the CM-3600d can be set to Soft-Flash mode. Numerical UV Control technology makes faster, more effective and reliable measurements of FWA treated materials such as textiles, papers and detergents.




Numerical Gloss Control

Numerical UV Control







The monolith polychromator unit assures highest durability and stability. With one exception (Observing Lens), the CM-3600d does not use any moving parts. Furthermore, each CM-3600d passes a hard endurance test program to comply with KONICA MINOLTA's highest quality standards

in accordance with ISO-9001

Performance

- Fast, simultaneous measurement of Specular Included and Specular Excluded components (SCI/SCE)
- Fast, Instantaneous numerical UV adjustment enables UV-included, UV-excluded, and UV-adjusted data to be obtained simultaneously
- Precise inter-instrument agreement.
 All CM-3600d's meet published inter-instrument agreement.


Versatility

- Full wavelength range 360-740nm with 10nm pitch
- Large (6") sphere; d/8° geometry
- Reflectance and transmittance sample measurements
- Changeable measurement areas (ø4mm,ø8mm,or ø25.4mm)
- View finder design for easy sample viewing
- Compact and lightweight

Reliability


High reliability design with fewest moving parts of any benchtop spectrophotometer


SpectraMagic™NX (Optional) Supports Windows®2000/XP



SpectraMagic™NX enables you to perform comprehensive color inspection and analysis of incoming raw materials, in process production, and outbound color critical goods and materials in virtually any industry. With SpectraMagic™NX you can insert digital images with measured data. Measure samples in any of 8 universally accepted color spaces. Select from 15 illuminants, and up to 40 indices to determine specific color and appearance properties, such as strength, brightness, haze, yellowness, opacity and strength. You can even configure up to 3 customized color equations. Reports range from simple Pass/Fail to trend charts. histograms, color plots, and spectral graphs. SpectraMagic™NX comes with predefined templates using skin technology, or you can create your own templates. For illustrations and explanations to understanding color and color measurement technology, there is a link to Konica Minolta's well known and respected "Precise Color Communication". Step by step navigation help. SpectraMagic™NX conforms to FDA 21 CFR Part 11 assuring integrity and reliability of data records

 Windows[®] is a trademark of Microsoft Corporation in the USA and other countries.





Specifications

Specifications	
Illumination/ observation system	Reflectance;d/8 (diffused illumination, 8-degree viewing), equipped with simultaneous measurement of SCI (specular component included) / SCE (specular component excluded) Conforms to CIE No.15,ISO7724/1,ASTME1164, DIN5033 Teil7 and JIS Z8722 condition C standard. Transmittance: d/0 (diffused illumination, 0-degree viewing) Conforms to CIE No.15, ASTME1164 and DIN5033 Teil7 standard.
Light-receiving element	Silicon photodiode array (dual 40 elements)
Spectral separation device	Diffraction grating
Wavelength range	360 to 740nm
Wavelength pitch	10nm
Half bandwidth	Approx.10nm
Reflectance range	0 to 200%; resolution: 0.01%
Sphere size	ø152mm
Light source	Pulsed xenon lamps (X4)
Measurement time	Approx. 1.5 seconds
Minimum interval between measurements	Approx. 4 seconds;when SCI/SCE measured
Measurement/ illumination area	LAV : ø25.4mm/ø30mm MAV: ø8mm/ø11mm (Selectable) SAV : ø4mm/ø7mm
Repeatability	Spectral reflectance: Standard deviation within 0.1% Colorimetric values: Standard deviation within ∆E*ab0.02
Inter instrument	Mean ΔE*ab0.15 (LAV/SCI) Average for 12 BCRA Series II color tiles
agreement	compared to values measured with master body.
Temperature dependence	Spectral reflectance: Within $\pm 0.10\%$ °C Color difference: Within ΔE^* ab 0.05/°C
UV adjustment	Instantaneous numerical adjustment
UV cut filter	400nm cutoff and 420nm cutoff
Transmittance chamber	Width: 133mm; depth: approx. 50mm; measurement dia.: approx. 17mm Transmission sample holder (Optional accessory): Sample holder; for both plate-shaped and liquid samples (removable)
Interface	RS-232C, D-SUB 9-pin (female) terminal
Power	AC120V/230V 50/60Hz (Using included AC adapter)
Operating temperature/ humidity range (*1)	13 to 33°C, relative humidity 80% or less (at 33°C) with no condensation
Storage temperature/ humidity range	0 to 40°C, relative humidity 80% or less (at 33°C) with no condensation
Size (WxHxD)/weight	244 x 208 x 378 mm(9-5/8 x 8-3/16 x 14-7/8 inch),12 kg(26-7/16 lb.)
Standard accessories	White Calibration Plate, Target Mask (ø4mm) , Target Mask (ø8mm) , Target Mask (ø25.4mm) , Zero Calibration box, AC adapter, Dust Cover, Accessory Case, Unit Driver, RS-232C Cable (9-pin,2m)
Options	SpectraMagic"NX(software) ,Transmittance Specimen Holder, Cell (2mm) / (10mm) / (20mm), Transmittance Zero Calibration Plate, RS-232C Cable (IBM,PC/AT 5m) / (IBM,PS/2 2m) / (IBM,PS/2 5m)

*1 Operating temperature/humidity range of products for North America: 13 to 33°C, relative humidity 80% or less (at 31°C) with no condensation

SAFETY PRECAUTIONS

For correct use and for your safety, be sure to read the instruction manual before using the instrument.




 Always connect the instrument to the specified power supply voltage Improper connection may cause a fire or electric shock.

Network construction for color control either within an organization or between organizations

High inter-instrument agreement between the same Konica Minolta model and also among all CM models (benchtops and portables): CM-2000 series, CM-3000 series,

This inter-instrument agreement is ideal when multiple units will be used for color control either within an organization or between organizations.



KONICA MINOLTA SPECTROPHOTOMETER LINEUP

CM-3700d series CM-3600d CM-3500d

Laboratory and Production modell Unique Top Port bench-top model CM-2600d/2500d/2500c Top class accuracy portables





"State of the Art" Reference models

Certificate No : YKA 0937154 Registration Date : March 3, 1995

Certificate No : JQA-E-80027 Registration Date: March 12, 1997

KONICA MINOLTA SENSING, INC. Konica Minolta Sensing Americas,Inc Konica Minolta Sensing Europe B.V.

New Jersey, U.S.A. European Headquarter /BENELUX German Office (International) German Office (Germany) French Office UK Office Italian Office Swiss Office Nordic Office Austrian Office

Osaka, Japan

Polish Office Konica Minolta (CHINA) Investment Ltd. SE Sales Division SE Beijing Office SE Guangzhou Office

Konica Minolta Sensing Singapore Pte Ltd.

KONICA MINOLTA SENSING, INC. Seoul Office **Phone:** 888-473-2656(in USA), 201-236-4300(outside USA)

Singapore

Seoul, Korea

Nieuwegein, Netherland Langenhagen, Germany München, Germany Roissy CDG, France Milton Keynes, United Kingdom Milan, Italy Dietikon, Switzerland Västra Frölunda, Sweden Wien Austria Phone: +31(0)30 248-1200 Phone: +49(0)511 7404-862 Phone: +49(0)89 630267-20 Phone: +33(0)1 493-82519 **Phone :** +44(0)1908 540-622 **Phone :** +39(0)23 90111 Phone: +41(0)43 322-9800 Phone: +46(0)31 7099464 Phone: +43(0)1 87882-430 Wien, Austria Warszawa, Poland Phone: +48(0)22 56033-00 Phone: +86-021-5489 0202 Phone: +86-010-8522 1551 Shanghai, China Guangzhou, China Phone: +86-020-3826 4220

Fax: +31(0)30 248-1211 Fax: +49(0)511 7404-807 Fax: +49(0)89 630267-67 Fax: +33(0)1 493-84771 Fax: +44(0)1908 540-629 Fax: +39(0)23 9011219

Fax: 201-785-2480

Fax: +41(0)43 322-9809 Fax: +46(0)31 474945 Fax: +43(0)1 87882-431 Fax: +48(0)22 56033-01 Fax: +86-021-5489 0005 Fax: +86-010-8522 1241 Fax: +86-020-3826 4223 Fax: +65 6560-9721

Phone: 02-523-9726 Fax: 02-523-9729 http://konicaminolta.com/about/se/contact.html

Phone: +65 6563-5533